

INVACARE® TOP END®
CROSSFIRE™
SERIES WHEELCHAIRS



The **Invacare® Top End® Crossfire™ Series** of everyday wheelchairs offers a range of frame materials and frame styles to suit a rider's individual needs and preferences. From the highly adjustable Invacare® Top End® Crossfire™ T7A Wheelchair to the Invacare® Top End® Crossfire™ T6 Wheelchair with a fixed rear seat-to-floor height, you can find the ideal chair for the most energy efficient ride. The Crossfire Series Wheelchairs are the perfect fit for anyone, from the new rider to the more experienced, who wants a custom, super light chair built for performance.

The **Invacare® Top End® Crossfire™ T7A Wheelchair** features an adjustable, performance-minded, open frame design unlike any other rigid chair available today.

Features

- Super light 7005 aluminum frame material
- Modular frame design allows for ample seat-to-floor height adjustments without re-squaring the casters
- 16 degrees of back angle adjustments
- 3.25" of center-of-gravity adjustments
- No charge mesh back upholstery offers more comfort in various climates
- Pressed bearing allows for a more controlled rotation and eliminates caster flutter

HCPCS Code: K0005



1. Adjustable angle fold down back
2. Invacare® Matrx® Seating available



Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Back Angle	Front Frame Angle	Camber	Tubing Diameter
CT7A	Seat Width + 6.75"	12"-20" W 14"-20" D	16"-21" Front 14"- 20" Rear .5" Increments	75°, 79°, 83°, 87°, 91°	70°, 75°, 80°, 85°	0°, 3°, 6°	1.12"
	Knee-to-Heel Range	Center-of-Gravity	Wheels	Casters	Total Weight**	Weight Capacity	Limited Warranty
	11.25"-20.25"	0.75"-3.75"	22", 24", 25", 26"	3", 4", 5"	17.8 lb.	250 lb.	5 Years Frame

**16"x16" with 24" lightweight performance spoke wheels, Primo Hi-Pressure tires, screw-on seat upholstery, mesh back upholstery with adjustable tension, long seat support and 4" urethane casters.

Please see the Owner's Manual or Invacare Configurator for complete option/specification combinations.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

The **Invacare® Top End® Crossfire™ T6 Wheelchair** offers a minimalist open frame design along with providing riders with the key adjustments such as center-of-gravity and back angle to keep the chair light and virtually maintenance free. In addition, custom options such as an Ergo seat or fixed back can be added to offer a more customized fit for the experienced rider.

Features

- Mono-tube, open frame design making vehicle transfer easy
- Choice of fold down or fixed back
- 4.25" of center-of-gravity adjustment
- Full offerings of wheel and tire options

HCPCS Code: K0005



1. Ergonomic seat option
2. Folding side guard
3. Mesh upholstery
4. Custom short fork
5. No caster caps

Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Back Angle	Front Frame Angle	Camber	Tubing Diameter
CRF	Seat Width + 6.75"	12"-20" W 14"-20" D	16"-21" Front 14"- 21" Rear .5" Increments	75°, 79°, 83° 87°, 91°	70°, 75°, 80° 85°, 90°	0°, 3°, 6°, 9°	1.12"
	Knee-to-Heel Range	Center-of-Gravity	Wheels	Casters	Total Weight*	Weight Capacity	Limited Warranty
	11.25"-20.25"	1.25"-5.5"	22", 24", 25", 26"	3", 4", 5"	19.5 lb.	250 lb.	5 Years Frame

* 16"x16" with 24" wheels, folding back, adjustable tension upholstery and 4" casters.

Please see the Owner's Manual or Invacare Configurator for complete option/specification combinations.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2014

Invacare Corporation.
All rights reserved.
Trademarks are
identified by the
symbols ™ and ®.
All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise
noted. Specifications
are subject to change
without notification.
Form No. 10-070
Rev. 01/14
131212



Yes, you can.®
