# WHEELCHAIRS, SCOOTERS, TRANSPORT CHAIRS & RAMPS

## **‡EZee Life**

#### PRODUCT DESCRIPTION, FEATURES & SPECIFICATIONS





### **OUR MISSION**

Since 1993, HPU Medical Wholesale has dedicated itself to developing and providing health care solutions that allow individuals, with physical challenges, the opportunity to safely enjoy the freedom of increased mobility and independence.

## **OUR VALUES**

We continually work closely with industry professionals to stay informed of market trends and our clients changing needs with an ongoing commitment to our four core values:

Continuously Improve—the quality of life for clients through increased comfort, function and mobility.

Be Accountable—for the solutions we provide.

Exceed Expectations—with results that match the needs of the individual.

Be Progressive—by developing, and introducing new products that reflect our level of dedication to improved mobility and independence.

## OUR GOAL

With compassion, commitment, and core values intact, we strive to remain a reliable and responsible source of better living solutions.

## CONTENTS

	Power wheelchairs
6	Wheelchairs.
1	Transport chairs.
1	Scooters.
22	Single Fold Ramp with Grip Tape.
23	Single Fold Wheelchair Ramps.
24	Multi-Fold Wheelchair Ramps.

#### **POWER WHEELCHAIRS**

#### CH4050 EZEE FOLD STANDARD

The EZee Fold Standard power wheelchair is a lightweight portable power wheelchair that folds easily for easy transport and storage but offers a range of up to 20 kms per battery. Optional second and third batteries will increase range accordingly.

With a seat width of 17 3/4" and a seat depth of 15 3/4" and a weight capacity of up to 264 lbs. this chair will fit most average sized adults.



6 kph (3 3/4 mph)

5 hours per battery

670 x 560 x 300 mm (26.4" x 22" x 11.8")

760 x 650 x 400mm (29.9" x 25.6" x 15 3/4")

15-20 km/battery

12 degrees

Max speed:

Endurance:

Overall size:

Folded size:

Max. Climbing Grade:

Charging period:

#### CH4060 EZEE FOLD HD

The EZee Fold HD power wheelchair is a lightweight portable power wheelchair that folds easily for easy transport and storage but offers a range of up to 10 kms per battery. Optional second and third batteries will increase range accordingly.

With a seat width of 17 3/4" and a seat depth of 15 3/4" and a weight capacity of up to 396 lbs. this chair will fit most larger sized adults.




450mm (17 3/4")
470mm (18 1/2")
400mm (15 3/4")
Height adjustable backrest.
25kgs (55 lbs) with 1 battery
180kgs (396 lbs)
8" x 2" Shock absorbing rubber caster
12" x 2" Wheels
Lithium battery, 24V, 7AH
6 kph (3 3/4 mph)
12 degrees
5 hours per battery
10 km/battery
918 x 603 x 892 mm (36.1" x 23.7" x 35.1")
668 x 603 x 324 mm (26.3" x 23.7" x 12.8")

#### WHEELCHAIRS

#### CH1090

The EZee Life™ Economy Wheelchair is available in either 16" or 18" seat width, has a 16" seat depth, an adjustable seat height, removable footrests and folds for easy storage and transportation. This wheelchair has a weight capacity of 300lbs. and is suitable for short trips and

occasional use.




35 1/2"

25 1/2"

Overall height: Overall width:

#### CH1091


This 22" heavy duty wheelchair is an economical choice for users who require a larger than standard wheelchair with a capacity of up to 400 lbs. The CH1091 features flat free tires, removable armrests and can be fitted with optional elevating leg rests. The rugged design of these wheelchairs makes them last for many years.



#### WHEELCHAIRS


#### CH1092

This EZee Life bariatric wheelchair is a heavy duty wheelchair with a seat width of 24", a weight capacity of 400lbs., steel folding frame, removable arm rests, removable foot rests, maintenance free solid tires and embossed light weight pvc upholstery. The rugged design of these wheelchairs makes them last for many years with only minor maintenance issues common among all wheelchairs to deal with.



#### CH1093

This EZee Life Standard Wheelchair is available with a 18" seat with and are an economical choice for those who need a wheelchair for regular to heavy use and institutional settings such as stores, hospitals, government offices etc. Because this standard wheelchair can be propelled by the user or pushed by a care giver, it will provide both independent and dependent mobility.



#### **WHEELCHAIRS**

#### CH1094


The EZee Life™ 20" standard wheelchair is a steel framed folding wheelchair with removable foot rests, removable arms, maintenance free solid tires, mag wheels and embossed light weight pvc upholstery. With their rugged design standard wheelchairs last for many years of service with only minor maintenance issues common among all wheelchairs to deal with.

SPECIFICATIONS		
Seat Width:	20"	
Seat Depth:	16"	
Seat Height:	17 1/2" or 19 1/2"	
Back height:	15 1/2"	
Arm height:	9 1/2"	
Arm length:	11"	
Foot rest length:	15" to 19"	
Footrest hanger angle:	60 degree	
Wheelbase: 18"		
Wheel size: 24"		
Caster size:	8"	
Tire type:	Solid	
Upholstery type:	PVC Vinyl	
Weight Capacity:	300 pounds	
Frame Material:	Steel	
Weight of chair:	37 lbs	
Overall height:	35 1/2"	
Overall width:	28"	
VICAL WHOLESALE		

#### CH1096

The 18" CH1096 Lightweight Wheelchair is constructed of an aluminum frame and has many standard features that other wheelchairs in this category offers as extra cost options. The dual cross brace frame style will provide a higher weight capacity and a structurally stronger frame than a single cross brace design.

This wheelchair's lightweight and conventional folding method will make it an easy chair to lift into vehicles for caregivers.




#### WHEELCHAIRS

#### CH1097

The 20" CH1097 Lightweight Wheelchair is constructed of an aluminum frame and has many standard features that other wheelchairs in this category offers as extra cost options. The dual cross brace frame style will provide a higher weight capacity and a structurally stronger frame than a single cross brace design.

This wheelchair's lightweight and conventional folding method will make it an easy chair to lift into vehicles for caregivers.



#### CH1041 - 17" TRANSPORT CHAIR

The EZee Life™ CH1041 transport chair is a lightweight wheeled chair that folds up for easy transportation and used for transporting people who are unable to walk distances required. While weighing only 19 pounds the CH1041has a weight capacity of 300 lbs. The Ch1041 features a seat width of 17", a fold down back, 8" flat free tires, and swing-away removable footrests with flip up footplates.



SPECIFICATIONS	
Seat Width:	17"
Seat Depth:	16"
Seat Height:	19"
Back height:	17"
Arm height:	10"
Arm length:	14 1/2"
Foot rest length:	8-9 1/2"
Caster size:	8"
Weight Capacity:	300 pounds
Frame Material:	Aluminium
Weight of chair:	19 lbs
Weight of footrests:	3.2 lbs
Overall width:	20"
Folded width:	10 1/2"

#### **TRANSPORTERS**

#### CH1042 - 19" TRANSPORT CHAIR

The EZee Life™ CH1041 transport chair is a lightweight wheeled chair that folds up for easy transportation and used for transporting people who are unable to walk distances required. While weighing only 20 pounds the CH1041has a weight capacity of 300 lbs. The Ch1041 features a seat width of 17", a fold down back, 8" flat free tires, and swing-away removable footrests with flip up footplates.



SPECIFICATIONS	
Seat Width:	19"
Seat Depth:	16"
Seat Height:	21"
Back height:	18"
Arm height:	10"
Arm length:	14 1/2"
Foot rest length:	8 1/2"
Caster size:	8"
Weight Capacity:	300 pounds
Frame Material:	Aluminium
Weight of chair:	20 lbs
Weight of footrests:	3.1 lbs
Overall width:	22"
Folded width:	9 1/2"

#### CH1044 – FEATHER LITE TRANSPORT CHAIR

The EZee Life™ CH1044 Featherlite transport chair features a high 325 pound weight capacity, flip back padded armrests and a lightweight aluminum frame that weighs 13.7 kg (30.2 lb) Its large 12" rear wheels and 8" casters makes this the easiest transport chair to push over rough terrains and obstacles while carrying larger adults. The standard breathable mesh back and seat pad allows the user's skin to breath and provide comfort when sitting for longer periods of time.





SPECIFICATIONS	
Seat Width:	18"
Seat Depth:	16"
Seat Height:	18"
Back height:	17 1/2"
Arm height:	10"
Arm length:	13"
Caster Size:	8"
Wheel Size:	12"
Weight Capacity:	325 lbs
Frame Material:	Aluminum frame
Weight of chair:	30.2 lbs
Overall height:	37 1/2"
Folded height:	29"
Overall width:	37 1/2"
Folded width:	11"
Overall depth:	37"
Folded depth:	41"

#### **TRANSPORTERS**

#### CH1043 - TRAVELER TRANSPORT

The Traveller is an 22.3 lb lightweight transport chair made of aluminum and steel frame. The Traveller will meet all your destination needs in airport terminals, shopping centres and many more locations. This lightweight model features a unique design that allows folding and carrying in its included carry bag to be quick and convenient.



SPECIFICATIONS	
Seat Width:	17"
Seat Depth:	16"
Seat Height:	19"
Back height:	16 1/2"
Arm height:	10 1/2"
Arm length:	10"
Caster Size:	8"
Wheel Size:	8"
Weight Capacity:	220 lbs
Frame Material:	Aluminum frame
Weight of chair:	22.3 lbs
Overall height:	37"
Folded height:	30"
Overall width:	22 1/4"
Folded width:	11 1/2"
Overall depth:	38"
Folded depth:	18"

#### **EZ5S SCOOTER**


The EZ5S mobility scooter is a 3 wheeled portable scooter small enough to be easily loaded into vehicles but powerful enough to be used both indoors and outdoors on solid and relatively smooth surfaces. The EZ5S Scooter disassembles into four pieces without the need of any tools for easy transport and storage with the rear section being the heaviest piece weighing only 31 lbs. (14kgs). With a range of up to 10 miles (16 kms) and a speed of up to 3.5 mph (5.6 kph) the EZ5S mobility scooter will last most users the entire day without the need of recharging. Range and speed will vary depending on user weight, battery condition and terrain.

SPECIFICATIONS	
Drive type:	Rear wheel drive
Maximum Speed:	3.5 mph
Maximum Climbing Slope:	8.00%
Range:	10 miles
Turning Radius:	31.5"
Ground Clearance:	2"
Suspension:	None
Weight capacity:	250 lbs
Battery:	12V 12 AH x 2, Sealed, deep cycle
Charging Time:	10 Hours
Rear Wheels:	7.5" x 2.5" Solid Tires
Front wheels:	7.5" x 2" Solid Tires
Overall Length:	41"
Overall Width:	19.5"
Folding:	no
Disassembly:	Yes
Basket:	Yes
Seat Width:	16.5"
Seat Depth:	14.5"
Seat Swivel:	Yes
Lights:	Yes
Battery gauge:	Yes
Horn:	Yes
Signals:	None
Rearview Mirror:	None
Overall weight:	90.2 lbs

#### **SCOOTERS**

#### **EZ6S SCOOTER**

The EZ6S mobility scooter is a 4 wheeled portable scooter small enough to be easily loaded into vehicles but powerful enough to be used both indoors and outdoors on solid and relatively smooth surfaces. The EZ6S Scooter disassembles into four pieces without the need of any tools for easy transport and storage with the rear section being the heaviest piece weighing only 31 lbs. (14kgs). With a range of up to 10 miles (16 kms) and a speed of up to 3.5 mph (5.6 kph) the EZ6S mobility scooter will last most users the entire day without the need of recharging. Range and speed will vary depending on user



#### **EZEE TRAVEL SCOOTER**


The Ezee Life EZee Travel Electric scooters are trike style mobility scooters that are super light foldable and easily transportable, and so are the best choice for electric scooters for on-the-move-riders! Folded and stored for the night, it will be ready when you are the next morning. Slide the Ezee Travel into the back of a van with virtually no lifting effort and unfold it again in no time. There are no connectors, latches or mechanisms to open when folding or unfolding. It's that easy!

SPECIFICATIONS	
Drive type:	Rear wheel drive
Maximum Speed:	3.7 mph
Maximum Climbing Slope:	6 degrees (10 %)
Range:	up to 6 miles
Turning Radius:	29"
Ground Clearance:	3.6"
Suspension:	No
Weight capacity:	250 lbs
Battery:	7 A/h 6 kg (13 lbs) or 12 A/h 8 kg (18 lbs)
Charging Time:	7 – 10 hours
Rear Wheels:	200 x 50 Solid tires (8" x 2")
Front wheels:	200 x 50 Solid tires (8" x 2")
Overall Length:	39"
Overall Width:	22"
Folding:	Yes
Disassembly:	Partial
Basket:	Underseat tray
Seat Width:	18"
Seat Depth:	16"
Seat Swivel:	No
Lights:	No
Battery gauge:	Yes
Horn:	No
Signals:	No
Rearview Mirror:	No
Overall weight:	72 lbs

#### **SCOOTERS**

#### EZEE CLASSIC SCOOTER

The Ezee Life EZee Classic Electric scooters are small lightweight mobility scooters that are super light, easily disassemble into four parts or can be folded for easy transportation and provide an excellent choice for electric scooters for on-the-move-riders! Because this mobility scooter easily breaks down into pieces it can be transported in almost any vehicle yet is strong enough to carry adults up to 250 pounds (115 kgs) in weight!



#### **F7FF ROVER 4 SCOOTER**

EZee Life EZee Rover 4 electric scooters are Rugged, Comfortable and Safe. Finally, an electric mobility scooter scooter that is not afraid of the outdoors and can take you places no other scooter can! The Ezee Rover 4 comes apart and easily folds for transport or storage and fits into the trunk of most cars. This Rugged mobility scooter comes with a full light package for optimal safety, and a keyed ignition for added security. A nice wide 18 inch deluxe memory foam seat, and big 10" pneumatic EASY-RIDE tires on the front wheels, and 14" tires on the rear wheels give the Rover 4 scooter a stable comfortable ride on most terrain



#### SINGLE FOLD RAMP WITH GRIP TAPE

#### 3' SINGLE FOLD RAMP W/ GRIP TAPE

The EZee Life™ single-fold EZ3GT wheelchair ramps provide indoor and outdoor accessibility for wheelchairs and scooters in a variety of environments. Built with the quality you expect, EZee Life™ ramps are lightweight but strong, available in a variety of length and feature grip tape for a non-slip surface.

- Lightweight
- Grip tape traction control
- Aluminum construction.
- Can be used indoors or outdoors
- Slip resistant
- Easy to install
- Maintenance-free
- 600 lbs Capacity







SPECIFICATIONS	
Width:	28.3"
Length:	3'
Maximum Rise:	6" with occupied device
	9" with unoccupied device
Construction:	Aluminum
Weight capacity:	600 lbs
Product weight:	16 lbs

#### SINGLE FOLD WHEELCHAIR RAMPS

#### **BI-FOLD WHEELCHAIR RAMP**

The EZee Life™ single-fold wheelchair Ramps provide indoor and outdoor accessibility for wheelchairs and scooters in a variety of environments. Built with the quality you expect, EZee Life™ ramps are lightweight but strong, available in a variety of lengths.

- Slip resistant
- Easy to install
- Maintenance-free
- 600 lbs Capacity
- Lightweight

- Textured surface traction control
- Aluminum construction
- Can be used indoors or outdoors





SPECIFICATIONS		
Width:	28.3"	
Length:	3', 4', 5', 6'	
Construction:	Aluminum	
Weight capacity:	600 lbs	
Product weight:	16 lbs, 22 lbs, 28 lbs, 34 lbs	



MAXIMUM RISES FOR BI-FOLD RAMPS		
Ramp Length	Occupied Device	Unoccupied Device
3'	6"	9"
4'	8"	12"
5'	10"	15"
6'	12"	18"

Note: Device refers to wheelchair or scooter.

#### MULTI-FOLD WHEELCHAIR RAMPS

#### MULTI-FOLD WHEELCHAIR RAMP

The EZee Life™ Multi-fold wheelchair ramps are portable ramps that provide access to areas up to 10" high while being approximately half of their extended width and length when folded. This folding design allows them to be stored away or transported in fairly small areas.

- Slip resistant
- Easy to install
- Maintenance-free
- 600 lbs Capacity
- Lightweight

- Textured surface traction control
- Aluminum construction.
- Can be used indoors or outdoors





SPECIFICATIONS		
Width:	29"	
Length:	5', 6', 7', 8', 10'	
Construction:	Aluminum	
Weight capacity:	600 lbs	
Product weight:	30 lbs, 35 lbs, 44lbs, 51 lbs, 63 lbs	



MAXIMUM RISES FOR MULTI-FOLD RAMPS		
Ramp Length	Occupied Device	Unoccupied Device
5'	10"	15"
6'	12"	18"
7'	14"	21"
8'	16"	24"
10'	20"	30"

Note: Device refers to wheelchair or scooter.

"EZee Life products are manufactured using modern manufacturing standards and equipment."



"All EZee Life products are factory quality tested to ensure reliability, soundness and safety."



