

**STELLATO**<sup>TM</sup>

# FOLDING WHEELCHAIR



**FUTURE**  
MOBILITY HEALTHCARE INC.

# STELLATO™

## Adjustable Axle and Seat Height

- Dual position axle plate with multiple vertical mounting holes
- Rear axle mount can be adjusted for greater stability and movement with seat height adjustment



## Adjustable Back Canes

- Adjustable height back canes can be set between -7 and 21 Degree



# STELLATO

## Flip Back Arms

- Flipback back arms to allow for easy transfer
- Choice of either fix for height adjustable arm



## Features on the Stellato Wheelchair

- Elliptical heavy duty crossbrace tubing
- Adjustment of caster housing angle
- Removable rear wheels with stainless steel or plastic coated handrims
- 60° or 70° swingaway upgraded front riggings
- Aluminum or composite foot plates
- Slip on seat upholstery
- Armrest options flipback or T-Style
- Wheel locks (push or pull to lock) with extension handles
- Adjustable rear anti-tips
- Quick release axles (optional)
- Heel loops
- Single aluminium cross-brace for optimum rigidity and ease of folding
- Elliptical heavy duty cross brace tubing creates a structurally sound frame

## Adjustable Anti Tippers

- Fully adjustable for various seat to floor height.



# ATTO

## SPECIFICATIONS:

SEAT WIDTHS	14"- 20", 22"
SEAT DEPTHS	15"- 20"
OVERALL SEAT WIDTH	Seat Width + 9"
SEAT TO FLOOR HEIGHT	Small Frame 13 1/4"- 18" Large Frame 15 1/2"- 19"
REAR WHEEL	20", 22" or 24" Composite Urethane
CASTER SIZE	4", 5", 6", 7", 8" Urethane
BACK STYLE	19" - 22" Low, 22" - 25" High Slip on back upholstery.
ADJUSTABLE BACK ANGLE	-7°, 0°, 7°, 14°, 21°
FRONT RIGGING	Pin Style, Swingaway, Dual Swingaway or Elevating Legrests
WEIGHT CAPACITY	250 lbs (114 Kg)
PRODUCT WEIGHT	31 lbs
FRAME MATERIAL	Lightweight Aluminum
ADP CODE	WMA2FM005

# CustomAnything!

*Customer Comes First.* Future Mobility Healthcare is a Canadian **Manufacturer** of Wheelchairs, Backs, Cushions and Accessories with our parts manufactured in house. This gives us the ability to customize anything needed for specific Customer Requests.



## WARRANTY:

The Stellato is warranted to be free from defects in materials and workmanship for a period of two (2) years from the date of purchase and lifetime for the frame. Upholstery for 90 days, upon normal usage by original purchaser.

Future Mobility Healthcare Inc. 3223 Orlando Drive, Mississauga, Ontario, L4V 1C5  
1-888-737-4011 Fax: 1-866-670-3377  
www.futuremobility.ca  
1 Buffalo River Place Buffalo, New York, 14210, USA  
1-855-9-FUTURE Fax: 716-783-9236,  
www.futuremobility.com  
All rights reserved

Printed in Canada.

Innovative Seating and Mobility Solutions



**FUTURE**  
MOBILITY HEALTHCARE INC.

