

UltraCare® XT

leading bed evolution







UltraCare® XT One Bed. A Total Solution.



It's a perfect bed for DELUXE ACCOMMODATIONS

The UltraCare XT has been created to look as great as it functions with features that offer residents the comforts of home. The UltraCare XT's UltraWide extension kit expands beyond a twin width and provides safer repositioning. For even more resident comfort, the PrevaMatt Elite surface is soft, yet highly resilient, gently cradling the resident's head and torso. Strategically-located mattress cuts promote easier flexing when the bed frame is contoured and keeps the resident in place to reduce the need for repositioning.



It's a superior bed for SKILLED NURSING FACILITIES

When gerontologist consultants guide design, the result is the best geriatric bed in the industry. Enhanced controls include features such as an intuitive pendant that is ergonomically functional with lighted buttons that enhance usability day or night. Add an assist handle with two positions for convenient use and perfect grip height to support resident transfer in and out of bed. We even partnered with Tempur-Pedic to offer mattresses that deliver enhanced pressure redistribution and comfort. When put all together, the end result is a versatile bed that meets your care challenges with safety in mind.



You Wanted Functional Speed and We Responded

QuickCare is standard, providing fast hi-lo adjustment.

Our QuickCare hi-lo speed raises and lowers the bed in less than 30 seconds. This increases bed operation efficiency under typical load applications and faster adjustments enable caretakers to deliver more enhanced resident care.

Better looking, better performing and more versatile, it is the only bed you'll ever need.



It's a flexible bed for BARIATRIC NEEDS

Flexibility is another design element where no other bed compares. The UltraCare XT can extend to a wider surface area, efficiently adding up to 7" to the sleep surface so you can more effectively manage your bariatric population. Also, with a safe working load of 500 lbs and a possible mattress width of 39" or 42", your Joerns bed offers more room for comfort and safer repositioning. AdvanceCare provides gravity-assist repositioning and comfort chair positioning, offering flexible care solutions while making repositioning safer and more comfortable for your residents and for your care providers.



It's an ideal bed for REHABILITATION NEEDS

The UltraCare XT is in a class of its own. The UltraLock option is designed for mobility and security at any height, promoting easier bed movement and safety. The four-section sleep surface has a full-range of comfort with therapeutic positions that improve circulation. QuickCare hi-lo speed promotes more efficient care while repositioning is made safer and more comfortable when using the 39" or 42" UltraWide extension kit. The trapeze can be added for rehabilitation needs and adjusts to accommodate most residents' reach, with three working positions for greater flexibility.



Unique Features to Help Reduce Resident Falls

77% of resident falls occur within the first 48 hours of admission. The **SafeLight** accessory provides soft floor lighting for residents to egress the bed during overnight hours.

The ergonomically mounted **Deluxe Assist Handles** provide a secure hand hold to support residents in maximizing their safety and comfort during bed transfers.

The evolutionary UltraCare XT from Joerns sets the standard. Once again.

